



Simplify. Streamline. Succeed.

0.9% Sodium Chloride Injection, USP

LESS WASTE, MORE EFFICIENCY:

With Assure Infusions means lower costs and a more reliable inventory



Products are preservative-free, PVC-free, DEHP-free, latex-free, and availability is subject to change. Please visit Assure Infusions products page at AssureIV.com for a current list of available services and new product announcements.

- ☒ Dependable Delivery: Offering fast, consistent, and high-quality product supply every time
- ☒ Clarity in Every Dose: Our advanced labeling features clear designs and prominent numbers to minimize medication errors and enhance safety
- ☒ Our cutting edge, cGMP compliant, and fully automated manufacturing process ensures a new benchmark for sterility, in a human touch-free environment
- ☒ Preservative-free, sulfite-free, latex-free, DEHP-free and PVC-free



Easy-to-read label with distinctive designs and prominent numbers to help reduce the potential for medication errors
cGMP-compliant compounded sterile products • Preservative-free, sulfite-free, latex-free, DEHP-free and PVC-free

Contact your local Assure Infusions Representative to learn how our products can benefit your facility.



863.292.2211



info@AssureIV.com

AVAILABLE PRESENTATIONS

Product	NDC#	Fill Volume	Container Type	Units/Case	Expiry Dating
0.9% Sodium Chloride Injection, USP	83155-100-25	250 mL	FFS	12	TBD
0.9% Sodium Chloride Injection, USP	83155-100-50	500 mL	FFS	12	TBD
0.9% Sodium Chloride Injection, USP	83155-100-10	1000 mL	FFS	12	TBD

Unmatched Automation

Our IV bags are manufactured, filled and terminally sterilized using advanced robotics in a human touch-free environment

The Pinnacle of Quality and Safety

Our pioneering approach eliminates human errors and contamination risks ensuring a precision level and consistency on every batch - far exceeding traditional 503b facilities

Advanced Visual Inspection

Every IV bag undergoes rigorous examination for particulates at a visual acuity surpassing the capability of human inspection

Unparalleled Capacity

With the ability to produce over two million bags per month, our facility is equipped to alleviate drug shortages and ensure a reliable, consistent supply

End to End Quality

- ☒ **Quality:** Products undergo rigorous testing and documentation processes, with final approvals by Quality in compliance with cGMP standards
- ☒ **Certificate of Analysis (COA):** Comprehensive test results are provided, including batch-specific COAs for our customers, demonstrating our commitment to surpass all 503B testing requirements
- ☒ **Sterility:** Assure has multiple methods of sterility testing, including ATP-bioluminescence when suitable, and membrane filtration, based on USP <1223> and USP <71> standards
- ☒ **Endotoxin Testing and Potency:** All batch lots are tested in strict accordance with USP <85> for Endotoxin and USP/NF monographs to ensure compliance with potency specifications
- ☒ **Environmental Monitoring:** Environmental monitoring is conducted throughout the compounding process, adhering to USP <1116> and <797> standards



Contact your local Assure Infusions Representative to learn how our products can benefit your facility.



863.292.2211



info@AssureIV.com